

### [Why are my protein marker bands not visible with the Lumitein™, One-step Lumitein™, or One-step Lumitein™ UV gel stains?](#)

Use of pre-stained protein markers can quench fluorescent protein stains. Switching to an unstained marker should resolve the issue.

### [Is protein staining on SDS-PAGE gels with the One-Step protein stains compatible with western blotting?](#)

We do not recommend staining proteins on the polyacrylamide gel with the One-Step protein gel stains before western blotting as it can significantly reduce protein transfer to the membrane. If protein detection on the gel prior to transfer is desired, the [Total Protein Prestains](#) would be more suitable.

The [VersaBlot™ Total Protein Normalization Kits](#) allow simple, sensitive and highly linear protein quantitation on SDS-PAGE gels and western blot membranes. The kits allow you to label purified proteins or cell lysates with our near-infrared [CF® dyes](#) before running the samples on SDS-PAGE. After electrophoresis, the bands can be visualized on the gel using a fluorescent gel scanner, allowing detection of as little as 1 ng protein per band. Labeled proteins also can be transferred to membranes for western blotting. The staining demonstrates excellent linearity for quantitation of total protein over a wide dynamic range, outperforming traditional western blot normalization based on housekeeping protein detection and staining is reversible.

### [Can the Lumitein™ Protein Gel Stain be used for staining native PAGE gels?](#)

Lumitein™ is a luminescent dye designed for detecting proteins in SDS polyacrylamide (SDS-PAGE) gels. It can also be used to detect proteins in native PAGE gels after an additional SDS incubation step as described in the protocol [PI-21002](#).